24 mesi-months-mois-monaten-meses

nome-name-nom-vorname-nombre

cognome-surname-prenom-nachname-apellido

indirizzo-address-adresse-auschrift-direccion

città-town-ville-ort-ciudad

SN _____ LP _____

data d'acquisto-purchase date-date d'achat
einkaufdatum-fecha de compra

versione-version-version-modell-versión

Timbro del rivenditore-Dealer's stamp-Cachet d'achat
Stempel der Fachhändler-Sello del revendedor

- FARO S.p.A. si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.
- FARO S.p.A. reserves the right to change the specifications of this equipment without notice.
- FARO S.p.A. se reserve le droit de modifier, sans préavis, les caractéristiques dans ce manuel.
- FARO S.p.A. behält sich recht vor, jederzeit stillschweigend technische oder bauliche Änderung vorzunehmen.
- FARO S.p.A. se reserva el derecho de modificar sin aviso previo las características incluidas en el presente manual de uso.

ED. 06/06 COD.988172 REV.1

FARO S.p.A. _____
Via Faro, 15
20060 Ornago (MI) Italia
Tel. : 039 / 6878.1 (r.a.)
Fax : 039 / 6010540

Internet: <http://www.faro.it>

MANUALE D'USO

POLILED COMPACT

MANUALE D'USO
INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRAUCHSANLEITUNG
MANUAL DE USO

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CE DECLARATION OF CONFORMITY

Il fabbricante FARO S.p.A. Via Faro n°15

The manufacture 20060 Ornago (Milano) ITALY

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Product's designation: Composite cure light

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It is newly manufactured in accordance with art. 6 par. 2 of Italian law D.L. n°. 626/94.

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It complies with the following European Directives and their implementations into national laws and subsequent modifications.

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CEI EN 60601-1 **CEI EN 60601-1-2**

It complies with the following technical standards:

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
GUARANTEE CERTIFICATE

FARO grants to the final customer a **24-months** guarantee starting from the date of the purchase invoice. Repairs under the guarantee shall be done at the **FARO premises or by an authorized FARO Service** person with transportation costs and risks at the expense of the purchaser. The guarantee certificate found on the endleaf of this booklet is considered valid only if it is entirely completed and bears the reseller's stamp. Repairs under guarantee (during the period of effectiveness) will be done only if the article to be repaired is accompanied by the shipping note or purchase invoice. The guarantee covers failures due to poor material quality or manufacturing defects and in case of legitimate claim the guarantee will allow repair or replacement free of charge. **No compensation for damages and/or interest will be recognized.** The guarantee is not deemed valid in the final judgement of FARO in case of tampering, damage, incorrect use, faulty maintenance or normal wear.

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SAFETY RULES

GB

Follow all the instructions set forth in this manual.

- This device is used to polymerize composite photohardening resins with photo starter activable in the band length between 440nm and 480nm in dentistry and must be used **in a dentist's office only by medical personnel**.
- Do not use in the presence of users of pacemakers or other implantable devices.
- Do not subject the gums or other soft tissues to the light beam. If necessary, protect any parts involved.
- Do not direct the light beam towards the eyes. Provide purposeful protection for very sensitive patients or those who undergo cataract surgery.
- Do not use this equipment for patients undergoing treatment with photosensitizing drugs. In case of risk, contact the specialized doctor.
- This device is not protected against water penetration.
- This device is not suited for use in environments with a mixture of anesthetic inflammable with air, oxygen or nitrous oxide.
- At the end of its operational life this device must be disposed of with separation of the different materials (plastic, metal, electronic circuits) following the legislation applicable for disposal of electronic equipment.
- Perform only the operations described in this manual and in all other cases apply to the Service Department.

FARO S.p.A. declines all responsibility for direct or indirect damage to persons or things in those cases where:

1. the instrument is not used as intended,
2. the instrument is not installed in accordance with the applicable standards and rules indicated,
3. the instrument is not used in accordance with all the instructions and rules set forth in this manual, and
4. personnel not authorized by FARO carry out modifications and repairs.

FEATURES

- The Poliled is used in dentistry as a polymerizer of photohardening composites.
- As a light source it uses a high-efficiency monochromatic blue LED.
- Wavelength emitted is between 440nm and 480nm.
- Although most composites are activated in this range of wavelengths, in case of uncertainty refer to the technical particulars of the compound or contact the producer.
- The electronic circuit **completely** integrated in the instrument is designed and built to ensure constant light emission.
- The Optical Fiber transmitting the light is easy to remove from the body of the device, can rotate 360° and is the **only component** which can and must be sterilized in an autoclave.

SAFETY INSTRUCTIONS DURING USE

- The Polymerizing Lamp is not supplied sterile, therefore it must be cleaned before use.
- Before using, always sterilize **the Optical Fiber only in the autoclave**.
- To avoid cross contamination between patients, insert the sterile Optical Fiber while taking care to avoid contact with contaminated hands or gloves.
- Do not direct the Polymerizing Lamp towards the patient's eyes.
- Turn off power before any maintenance work.
- Before any treatment, make sure the instrument works perfectly.
- Always insert the Optical Fiber completely to avoid weakened yield.
- Replace damaged Optical Fibers without delay.

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INSTALLATION**Make sure the following components are in the package.**

⇒ Poliled Cord.

⇒ Installation Manual.

SPECIALIZED TECHNICIANS MUST INSTALL THIS DEVICE. POWER MUST BE DISCONNECTED WHEN INSTALLING.

- Fasten the Power Cord inside the Chair using the Tearproof Cord so that a pull on the Cord will not cause traction on the power cord.
- Connect the Power cords only to a very low-voltage safety transformer (see Safety Instructions).

OPERATING INSTRUCTIONS

The Poliled Polymerizing Lamp has two different types of exposition, as follows.

- **Cycle Power:** 20-second exposition cycle at full power.
- **Soft Ramp Cycle:** 20-second exposition cycle with gradual increase in light intensity for first 10 second and up to full power.

USING AT FULL CYCLE POWER

- Press button (3) shortly. LED goes on full power for 20 seconds.
- A short acoustic signal is emitted at the beginning and halfway thru cycle. A longer sound is heard at cycle end.

USING WITH SOFT-RAMP FUNCTION

- Hold down button (3) for about 2 seconds until cycle start. The LED lights up and increases intensity gradually.
- A short acoustic signal is emitted at the beginning and halfway thru cycle. A longer sound is heard at cycle end.

VOLUNTARY STOPPING THE CYCLE

- Exposure can be stopped at any moment by again pressing a button (3) whether in Power or Soft-Ramp mode.

NOTE

- To obtain full output, make sure the Optical Fiber is inserted correctly and is not fouled with polymerizing material.
- In accordance with std 60601-1, the instrument is suited for intermittent use with 80 seconds ON and 10 minutes OFF.
- For extremely heavy applications (**approximately 15 consecutive cycles in Power mode**) a thermal safety protection shuts off the instrument. In these cases, any increase in temperature outside of the instrument is considered normal. Use will be inhibited for a few minutes. Complete resetting will take place after approximately 15 minutes.

MAINTENANCE**THE INSTRUMENT IS NOT PROTECTED AGAINST THE PENETRATION OF LIQUIDS. STERILIZE ONLY THE OPTICAL FIBER.****To sterilize the Optical Fiber, proceed as follows.**

- Withdraw the Optical Fiber from the body of the Polymerizing Lamp by pulling it axially.
- Sterilize the Optical Fiber in an autoclave.
- Refit the Optical Fiber by rotating it slightly down to the striker.

CLEANING**To clean the Lamp proceed as follows.**

- Turn off the Main Power Switch before cleaning.
- Use a soft cloth slightly moistened with alcohol or similar detergent while being careful not to allow alcohol to run inside the equipment.

- Any residual reconstruction material encrusted on the Optical Fiber or the Lamp Body must be removed with beveled plastic tools.
- To get highest output, always clean both sides of the Optical Fiber with a soft cloth every day. Both sides of the Optical Fiber must be cleaned before each sterilization.
- Periodically examine the condition of the Optical Fiber. Its surface should appear smooth, uniform and free of halos or spots. Any irregularities have a detrimental effect on polymerization.

TROUBLESHOOTING

The Polymerizing Lamp does not light.

- Make sure electric power cord is correctly plugged in.
- If the problem continues, call on our Service Department.

Polymerization power is slight.

- Check condition and cleanliness of Optical Fiber and change if necessary.
- If the problem continues, call on our Service Department.

Optical Fiber is lost or broken.

- To find, apply only to FARO Service centers.

Polymerizing Lamp shuts off too soon, emits three beeps and stops further cycles.

- The Thermal Protection has operated. Wait a few minutes before using again.

The surface of the instrument has heated up.

- You have performed many cycles. Let the instrument rest.
- If the problem continues, disconnect the instrument from the cord and cool it, then reconnect it. If the instrument heats up again without being used, completely disconnect it and apply to Service Department.

TECHNICAL PARTICULARS

Application: Polymerization or photohardening of dental materials

Input voltage, Vac: 17-26 /19-34 50 HZ

Consumption, VA: 11.5

Protection level against direct and indirect contact: device type B 

Wavelength: 440nm to 480nm

Protection against electric danger: device class II

Dimensions of instrument: Length 165mm

Cable: fixed with the instrument

Optical Fiber: dia. 8, autoclavable at 134°C Maximum dia.: 27mm

Applicable standards: EN 60601-1 and 60601-1-2

Weight: 110g

Light source: One 5W LED with optical concentrator

Working time:

Standard cycle at full power, seconds: 20

Ramp cycle from zero to full power, seconds: 10 - Ramp power at full power, seconds: 10

Acoustic Signals: Start and after every 10 seconds

This equipment while in its transportation storage packing can withstand the following environmental conditions for a period of 15 weeks at most:

⇒ **Surrounding temperature from -20°C to +70°C**

⇒ **Relative humidity from 10% to 90%**

⇒ **Atmospheric pressure from 500 to 1060 mbar.**

Disposal container. Dispose the device observing the standards about the differentiated collection of electrical material.

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